



Pathology and Laboratory Medicine



Orderable - GLUR GLUA (fasting glucose)

Turnaround Time: 4 hours STAT: 1 hour

Alternate Name(s):

Random Glucose Fasting Glucose

Specimen:



Laboratory: Core Lab



Requisition:GENERAL LABORATORY
REQUISITION



Method of Analysis: Enzymatic rate





Test Schedule: As required

Adult
Pediatric

4.5 mL Green top
Vacutainer tube
or
2 mL Grey top
Vacutainer tube
Vacutainer tube

2 mL Grey top
Vacutainer tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Minimum volume required for a pediatric sample is 0.2 mL

Reference Ranges:

Random	
<u>Age</u>	<u>Range</u>
< 1 year	2.6-11.0 mmol/L
> 1 year	3.4-11.0 mmol/L
Fasting levels - no caloric intake for at least 8 hours	
3 years - 16 years:	2.8-6.1 mmol/L
Adult (greater than 16 years):	3.9-6.1 mmol/L







Pathology and Laboratory Medicine

Critical Value	
<u>Age</u>	<u>Range</u>
0-3 days:	<2 mmol/L or >15 mmol/L
>3 days:	<2.8 mmol/L or >30 mmol/L

Interpretive Comments:

2018 CDA Updated Diagnostic Criteria for Diabetes:

FPG ≥7.0 mmol/L

Fasting = no caloric intake for at least 8 hours

0

A1C ≥6.5% (in adults)

0

2hPG in a 75 g OGTT ≥11.1 mmol/L

or

Random PG ≥11.1 mmol/L

Random = any time of the day, without regard to the interval since the last meal

Please check the following link for more details: http://guidelines.diabetes.ca/cpg

Special Processing:

Li Heparin (green top tube) unspun: stable for 2 hours

Li Heparin (green top tube) **spun**: stable for 8 hours at room temperature and 24

hours refrigerated.

Na Fluoride (grey top): spun and unspun: stable for 48 hours.