

Orderable - OXALU

Turnaround Time: 10 days

Specimen:

Random urine	
(24-hour urine is the preferred sample type; see Lab Test	
Information Guide entry for Citrate, 24-hour Urine)	

Collection Information:

Collect random urine in a sterile container and deliver to the core lab as soon as possible so that acid preservative can be added.

London Health Sciences Centre

Pathology and Laboratory Medicine

Minimum volume of urine required is 5 mL.

Patient should restrict intake of vitamin C-rich foods and vitamin C supplements for 48 hours prior to urine collection.

Reference Ranges:

Age	Range
< 6 months:	≤ 0.37 mmol/mmol
	creatinine
6-24 months:	≤ 0.26 mmol/mmol
	creatinine
2-5 years:	≤ 0.14 mmol/mmol
	creatinine
6-12 years:	≤ 0.08 mmol/mmol
	creatinine
> 18 years:	≤ 0.04 mmol/mmol
	creatinine

Interpretive Comments:

High values are found in primary hyperoxaluria, diseases of the distal ileum, pyridoxine deficiency, and ethylene glycol poisoning. Consumption of oxalate-rich foods (e.g. rhubarb, strawberries, spinach, and tomatoes) may increase excretion.



Requisition: GENERAL LABORATORY REQUISITION

Laboratory: Core Lab

Method of Analysis: Enzymatic

0_0

Test Schedule: Referred out Monday-Thursday



Referred Out Location: In-Common Laboratories



Pathology and Laboratory Medicine

OXALATE, RANDOM URINE

Special Processing:

Add 4 drops of 6 N HCl to random urine collection. Final pH should be 1-3.

Storage and Shipment:

Store and ship frozen.