

# MYOSITIS ANTIBODIES PROFILE

## Orderable - .MYOAP

Turnaround Time: 5 days

### Alternate Name(s):

Autoimmune Myositis Panel  
Myositis-specific antibodies  
Anti-Mi-2 alpha  
Anti-Mi-2 beta  
Anti-TIF1 gamma  
Anti-MDA5  
Anti-NXP2  
Anti-SAE1  
Anti-Ku  
Anti-PM-Scl 100  
Anti-PM-Scl 75  
Anti-JO1  
Anti-PL-7  
Anti-PL-12  
Anti-EJ  
Anti-OJ  
Anti-Ro-52  
Anti-Ha  
Anti-Ks  
Anti-Zo Anti-cN1A  
Anti-SRP



**Laboratory:**  
Clinical Immunology Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Immunoblot



**Test Schedule:**  
Batch analysis

### Specimen:


Adult	Pediatric
5 mL Gold top Vacutainer tube or 6 mL Red top Vacutainer tube	0-2 years: Red 0.5 Microtainer 2-10 years: 2 mL Red top
<b>EDTA, Heparin or Citrate plasma are also acceptable</b>	

### Collection Information:

Minimum volume of serum or plasma required for testing is 0.5 mL.

### Reference Ranges:

Negative



# MYOSITIS ANTIBODIES PROFILE

## Interpretive Comments:

Immunoblot results will be reported as either negative or positive for specific antibodies as the following.

- 1+ Weak Positive
- 2+ Positive
- 3+ Strong Positive

## Comments:

**This test is available exclusively to Rheumatologists and Respiriologists at LHSC/SJH and accepted from referred in locations.**

## Storage and Shipment:

Sample can be stored for up to 14 days at 4°C. For longer storage, store at -20°C.