

ANTIMONY,
URINE

Orderable - SBU / SBU24

Turnaround Time: 10 days

Specimen:

24-hour urine collected in an unused 24-hour urine container or random urine

Collection Information:

[Trace Elements Collection Requirements](#)

Reference Ranges:

SI Units

Random Urine:
0-108 nmol/mol creatinine

24 Hour Urine:
0-0.69 nmol/d

Conventional Units

Random Urine:
0-116 ng/g creatinine

24 Hour Urine:
0-0.084 µg/d

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information [here:](#)

Specimen Stability:

Ambient: 14 days
Refrigerated: 11 months
Frozen: 11 months



Laboratory:
Trace Elements Lab



Requisition:
TRACE ELEMENTS
REQUISITION



Method of Analysis:
ICP-MS/MS



Test Schedule:
Batched analysis