

ANTI-  
THYROGLOBULIN  
ANTIBODIES,  
SERUM/  
PLASMA

## Orderable - ATG

Turnaround Time: 24 hours

### Alternate Name(s):

Thyroid Antibodies

### Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 0.5 mL Red or Gold top Microtainer 2-10 years: 3 mL Red top Vacutainer tube
Lavender top (EDTA) tubes are also acceptable	
Light Green top (Li-Heparin) tubes are <b>not</b> acceptable	

### Collection Information:

Minimum volume of serum required is 700 µL for adult samples or 200 µL for pediatric samples.

### Reference Ranges:

≤ 115 IU/mL

### Interpretive Comments:

The presence of anti-thyroglobulin antibodies may cause falsely low thyroglobulin measurements.

### Comments:

Measurement of anti-thyroglobulin antibodies is automatically performed with each serum thyroglobulin test order.



**Laboratory:**  
Core Lab




**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Roche  
Electrochemiluminescence



**Test Schedule:**  
As required



ANTI-  
THYROGLOBULIN  
ANTIBODIES,  
SERUM/  
PLASMA

Thyroglobulin concentrations > 2000 µg/L may lead to falsely elevated anti-thyroglobulin antibody results.

**Storage and Shipment:**

Store and ship serum or plasma frozen.