

GABAPENTIN,
SERUM/
PLASMA

Orderable - GABA

Turnaround Time: 5 days

Alternate Name(s):

Neurontin

Specimen:

6 mL Red or 6 mL Green or 4 mL Lavender top Vacutainer tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Minimum volume required for testing is 0.5 mL.

Reference Ranges:

12.0 - 120.0 µmol/L

Interpretive Comments:

Alert value is 140.0 µmol/L.

Minimum effective concentration is 12.0 µmol/L.

A desirable clinical outcome may be expected between 70.0-120.0 µmol/L. Toxicity may occur with increasing frequency above 146.0 µmol/L

Storage and Shipment:

Separate serum and freeze. Ship sample to the lab cold or frozen.



Laboratory:
Toxicology/Special
Chemistry Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS



Test Schedule:
Monday--Friday
0800-1600