

HISTOPLASMA ANTIGEN

Orderable – HISTOAG

Turnaround Time: 2 days from receipt at reference laboratory

Specimen and Collection Information:

Please review the MiraVista Diagnostics Laboratory Test Information Index:

[Histoplasma Quantitative EIA Test](#)

Specimen Type	Minimum Volume	Collection Kit
Urine	0.5 mL	Sterile Container
Serum	1.2 mL	Gold Vacutainer (SST)



Laboratory:
Microbiology



Requisition:



Method of Analysis:
Quantitative Sandwich
Enzyme Immunoassay
(EIA)



Test Schedule:
Referred out weekdays



Referred to:
MiraVista Diagnostics
4705 Decatur Blvd.
Indianapolis, Indiana
46241
USA

Limitations

- Result is not intended to be used as the sole means for clinical diagnosis or patient management decisions
- Cross-reactions are seen with blastomycosis, paracoccidioidomycosis, penicilliosis, less frequently in coccidioidomycosis, rarely in aspergillosis and possibly sporotrichosis

Comments:

Testing restricted to Infectious Diseases (ID) service following a Microbiologist consultation.