

# **Orderable - CARB**

Turnaround Time: 24 hours

## Alternate Name(s):

Adult

Anticonvulsant Tegretol

# Specimen:



Laboratory: Core Lab

Requisition: GENERAL LABORATORY REQUISITION



Method of Analysis: Kinetic Interaction of Microparticles in Solution (KIMS)

#### 6 mL Red top Vacutainer tube 0-2 years: 0.5 mL Red top Microtainer 2-10 years: 4 mL Red top Vacutainer tube Green and gold top Vacutainers with separator gels are not acceptable.

Pediatric

**Collection Information:** 

Collect blood aseptically in a Vacutainer tube.



Test Schedule: As required

Other times by prior arrangement.

### **Reference Ranges:**

17 - 51 μmol/L

Critical Value: Toxic: >64 µmol/L

## Comments:

CRITICAL VALUE to be phoned to Nurse or Physician immediately.



Pathology and Laboratory Medicine