GROUP AND SCREEN

Orderable - GROUPANDSCREEN (LHSC/SJH)

Turn Around Time: 2 hours

STAT: 45 minutes

London Health Sciences Centre

Pathology and Laboratory Medicine

Specimen:

Age	Specimen Required
Pediatric <4 months:	1 mL Lavender (EDTA)
Infants >4 months and <3 years:	2 mL Lavender (EDTA) top
Children >3 years and	6 mL Pink (EDTA) top Vacutainer
Adults:	tube



Laboratory: Blood Transfusion Laboratory



Requisition: BLOOD TRANSFUSION REQUISITION or Electronic order

Reference Ranges:

See report



Method of Analysis: Anti-Human Globulin test

Test Schedule:

Dailv

Interpretive Comments:

(Used for Blood Grouping, Antibody Screen and Crossmatching)



<u>Comments:</u>

Crossmatched blood can be available within 5 to 10 minutes upon receipt of request to convert Group and Reserve to crossmatch.

Storage and Shipment:

Transport to the Blood Transfusion Laboratory within 24 hours of collection at room temperature.