





Pathology and Laboratory Medicine

Orderable - HpAlgGAb

Turn Around Time: 48 hours

Alternate Name(s):

Anti HAV IgG Hepatitis Anti HAV IgG

Specimen:

Adult	Pediatric
5 mL Gold top	0-2 years: 0.5 mL Gold top
Vacutainer	Microtainer
	2-10 years: 2 mL Gold top tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

See report

Interpretive Comments:

Indicates immunity either from past infection or immunization.

Comments:

Heparinzed samples should be avoided as they cause erroneous results.

0

Laboratory: Core Lab



Requisition:

GENERAL LABORATORY REQUISITION



Method of Analysis:

Electrochemiluminescence Immunoassay



Test Schedule:

24/7