

# HLA PRE SCREEN

## Orderable – HLA PRESCREEN-RENAL

Turnaround Time: 14 days

### Specimen:

2 x 4 mL Lavender (EDTA) top, 1 x 7mL Red top Vacutainer tubes

### Collection Information:

Samples that are being shipped into the laboratory can be drawn Monday- Friday

### Reference Ranges:

See report

### Special Processing:

Remove serum from clot. Keep cold.

### Storage and Shipment:

Keep EDTA blood at room temperature and deliver within 72 hours.



**Laboratory:**  
Transplant Lab



**Requisition:**  
[TRANSPLANT  
LABORATORY  
REQUISITION](#)



**Method of Analysis:**  
HLA Class I and Class II  
 typings are performed  
 using molecular  
 techniques (RSSO or SSP)



**Test Schedule:**  
As required Monday to  
 Friday except Stat  
 holidays