

Pathology and Laboratory Medicine

LAMOTRIGINE, SERUM/ PLASMA

# **Orderable - LAMICTAL**

Turnaround Time: 10 days

### Alternate Name(s):

Lamictal

#### Specimen:

6 mL Red top Vacutainer tube or 4.5 mL Green (Lithium Heparin) top Vacutainer tube

Collect blood aseptically in a Vacutainer tube.

Requisition: GENERAL LABORATORY REQUISITION

İ

Method of Analysis: LC-MS/MS

Laboratory:

Chemistry Lab

Toxicology/Special

**Reference Ranges:** 

**Collection Information:** 

4.0-39.0 µmol/L



**Test Schedule:** Once per week

## **Interpretive Comments:**

Risk of toxicity increases above 60 µmol/L

#### **Storage and Shipment:**

Separate and freeze. Ship on ice packs or frozen.