

SODIUM, URINE - RANDOM

Orderable - NAU

Turnaround Time: 4 hours

STAT: 1 hour

Specimen:

Random urine

Collection Information:

Collect a minimum of 5 mL of random urine in a sterile container.

Reference Ranges:

40-200 mmol/L

Interpretive Comments:

Increased sodium excretion may be found in adrenal cortical insufficiency, renal tubular disease and in patients receiving natriuretic agents.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Potentiometric, Indirect
Ion-Selective Electrodes



Test Schedule:
As required