

# TETANUS SEROLOGY

## Orderable - TETAN

Turnaround Time: Up to 28 days

### Specimen Type and Collection Information:

Please review the Public Health Ontario Laboratory Test Information Index:

[Tetanus Serology](#)

Specimen Type	Minimum Volume	Collection Kit
Whole blood or serum	5.0 mL blood or 1.0 mL serum	Blood, clotted – Gold vacutainer tube (SST)

Testing for Tetanus is not done for immunity/acute diagnosis/pre-vaccination screening.

### Storage and Shipment:

The specimen must be refrigerated.  
The specimen must be received within 24 hours of collection.



**Laboratory:**  
Microbiology



**Requisition:**  
See electronic order



**Method of Analysis:**  
EIA



**Test Schedule:**  
Referred out weekdays



**Referred to:**  
Public Health Laboratory  
1200 Commissioners Rd E  
London, Ontario  
519-455-9310