TOPIRAMATE SERUM/ PLASMA

Orderable - TOPI

Turnaround Time: 5 days

Alternate Name(s):

Торатах

Specimen:

Adult	Pediatric
6 mL Red or 6 mL Green	0-2 years: Red 0.5pk. or 0.6
or 4 mL Lavender top	Green Microtainer
Vacutainer tube	2-10 years: 2 mL Red or Dark
	Green top

Collection Information:

Collect blood aseptically in a Vacutainer tube with no separator.

Minimum volume required for testing is 1 mL of plasma or serum.

Reference Ranges:

15-60 µmol/L

Interpretive Comments:

Risk of toxicity increases above 75 μ mol/L.

Storage and Shipment:

Separate and freeze. Ship sample to the lab cold or frozen.



Laboratory: Toxicology/Special Chemistry Lab

GENERAL LABORATORY



REQUISITION

Requisition:

Å

LC-MS/MS

Method of Analysis:



Test Schedule: Monday--Friday 0800-1600



Pathology and Laboratory Medicine