



TRANSFERRIN
ISOELECTRIC
FOCUSING,
SERUM

Orderable - TIF

Turnaround Time: 3-4 weeks

Specimen:

5 mL Gold top Vacutainer is preferable

Collection Information:

All Transferrin Isoelectric Focusing testing requires approval to be performed. For approval, please contact Dr. Tony Rupar at 519-685-8500 ext. 71558 or pager 15880.

Minimum volume required is 1 mL.

Reference Ranges:

A descriptive report will be sent

Interpretive Comments:

Isoelectric focussing of transferrin is used as a test to detect Carbohydrate Deficient Glycoprotein Syndrome.

Storage and Shipment:

Send sample frozen on dry ice.



Laboratory:
Core Lab



Requisition:
[SICKKIDS METABOLIC
DISEASES LABORATORY
REQUISITION](#)



Test Schedule:
Referred out Thursdays



Referred Out Location:
Metabolic Diseases
Laboratory at Sick Kids
555 University Avenue,
Room 3333, Black,
Toronto, ON M5G 1X8,
Telephone: 416-813-7200