

VORICONAZOLE

Orderable - VOR

Turn Around Time: 10 days

Specimen:

Adult	Pediatric
4.5 mL Green (Lithium Heparin) top Vacutainer tube	2 mL Green top tube



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)

Collection Information:

Collect trough specimen immediately prior to next dose or at least 12 hours post-dose.

Minimum volume of plasma required for testing is 0.5 mL

Reference Ranges:

Therapeutic:	1.0 - 5.0 mg/L
Toxic:	> 6.0 mg/L(Trough)

Specimen Stability:

Refrigerated: 1 day

Frozen: >1 day

Storage and Shipment:

Store and send frozen.