

LACTATE CSF

Orderable - LACCSF

Turn Around Time: 4 hours

STAT: 1 hour

Specimen:

CSF

Collection Information:

Collect a minimum of 500 uL of CSF in a standard CSF tube.

CSF must **not** be collected into any anticoagulant container.

Keep sample cold when sending to laboratory.

Reference Ranges:

See report

Comments:

RBC hemolysis may cause interference (xanthochromia) with results.

Storage and Shipment:

Once collected, store and ship sample on ice to laboratory.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Endpoint methodology



Test Schedule:
As required