

# LACTATE, CSF

## Orderable - LACCSF

Turnaround Time: 4 hours

STAT: 1 hour

### Specimen:

CSF

### Collection Information:

Collect a minimum of 500 uL of CSF in a standard CSF tube.

CSF must **not** be collected into any anticoagulant container.

### Reference Ranges:

See report

### Comments:

RBC hemolysis may cause interference (xanthochromia) with results.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Endpoint methodology



**Test Schedule:**  
As required