

# LIDOCAINE SERUM/ PLASMA

## Orderable - LIDO

Turn Around Time: 1 week

### Alternate Name:

Xylocaine

### Specimen:

Adult	Pediatric
6 mL Red top Vacutainer tube	0-2 years: 0.5 mL Red top Microtainer 2-10 years: 3 mL Red top Vacutainer tube
Plasma (heparin or EDTA) is also acceptable	

### Collection Information:

Blood collection tubes with separator gels are not acceptable for this test.

Collect trough specimen within 1 hour prior to next dose.

Minimum volume of serum (or plasma) required is 200 µL.

### Reference Ranges:

1.5-5.0 mg/L

### Storage and Shipment:

Due to limited stability specimen must be received at ICL within 3 days of collection. Store and send cold.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Immunoassay



**Test Schedule:**  
Referred out Monday-  
Thursday



**Referred Out Location:**  
[In-Common Laboratories](#)