

MAPROTILINE,
SERUM/
PLASMA

Orderable - MAP

Turn Around Time: 5 days

Alternate Name(s):

Ludiomil

Specimen:

2 x 5 mL Gold top Vacutainer tube or 2 x 4 mL Lavender EDTA Vacutainer tubes

Collection Information:

Minimum volume of serum or EDTA plasma required is 3 mL.

To monitor therapy, draw trough specimen prior to next dose or 10 - 12 h after last drug administration. Separate as soon as possible.

Reference Ranges:

180 -720 nmol/L

Toxic: Greater than 2900 nmol/L

Storage and Shipment:

Separate serum as soon as possible.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)