

METHADONE,
SERUM

Orderable - METHAD

Turn Around Time: 14 days

Specimen:

Adult
6 mL Red top Vacutainer tube
6 mL Royal Blue top Vacutainer tubes for serum collection are also acceptable.

Collection Information:

Blood collection tubes with separator gels are not recommended for this test.
Minimum volume of serum required is 1.2 mL.
Biochemist approval is required for this test.

Reference Ranges:

Reference Range
100 – 400 ng/mL
Potentially toxic range: ≥ 2000 ng/mL

Storage and Shipment:

Store and ship serum frozen.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
GC/MS



Test Schedule:
Referred out Monday –
Thursday



Referred Out Location:
[In-Common Laboratories](#)