

 PAROXETINE,
SERUM

Orderable - PAROX

Turn Around Time: 7 days

Alternate Name(s):

Paxil

Specimen:

6 mL Red or Navy top Vacutainer tube

Collection Information:

Collect blood aseptically in a Vacutainer tube 10 - 14 hours after last dose.

Collect trough specimen prior to first morning dose.

A minimum volume of 0.5 mL of serum is required for testing.

Reference Ranges:

Up to 150 nmol/L

Storage and Shipment:

Separate, store and send frozen.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
HPLC



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)