

PIGEON IGG ANTIBODIES

Orderable - AAPG

Turn Around Time: 10 days

Alternate Name(s):

Pigeon Precipitins

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube or 6 mL Red top Vacutainer tube	0-2 years: 0.5 mL Red Microtainer 2-10 years: 2 mL Red top tube

Collection Information:

A minimum volume of 1 mL or serum is required for testing.

Grossly hemolyzed or lipemic samples not accepted.

Reference Ranges:

<30 mg/L

Interpretive Comments:

This assay detects antibodies to Pigeon serum proteins, feathers and droppings.

Specimen Stability:

Refrigerated: 7 days

Frozen: 3 months

Storage and Shipment:



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Fluorescence Enzyme
Immunoassay (FEIA)



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)



PIGEON IGG
ANTIBODIES

Store and send cold.