

PROPAFENONE  
SERUM/  
PLASMA

## Orderable - PROPA

Turn Around Time: 2 weeks

### Alternate Name(s):

Rythmol

### Specimen:

Adult	Pediatric
6 mL Red top Vacutainer tube	0-2 years: 2 x 0.5 mL Red top Microtainers 2-10 years: 3 mL Red top Vacutainer tube
Plasma (Sodium Heparin) is also acceptable.	

### Collection Information:

**Blood collection tubes with separator gels are not acceptable for this test.**

Collect specimen 2-6 hours post oral dose.

Minimum volume of serum (or plasma) required is 500 µL.

### Reference Ranges:

0.2-1.6 µg/mL

### Storage and Shipment:

Store and ship serum (or plasma) frozen.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Chromatography/mass  
spectrometry



**Test Schedule:**  
Referred out Monday-  
Thursday



**Referred Out Location:**  
[In-Common Laboratories](#)