

SODIUM, FLUID

Orderable - NAF

Turn Around Time: 4 hours

STAT: 1 hour

Specimen:

Fluid

Collection Information:

Collect a minimum of 5 mL of fluid in a sterile container.

Reference Ranges:

See report

Interpretive Comments:

Used in diagnosis and treatment of Aldosteronism, diabetes insipidus, adrenal hypertension, Addison's disease, dehydration, inappropriate ADH secretion, electrolyte imbalance.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Indirect Ion-Selective
Electrodes



Test Schedule:
As required