

THIOPURINE
METHYLTRANSFERASE:
GENOTYPE

Orderable - TPMTG

Turn Around Time: 21 days

Alternate Name(s):

TPMT Genotype

Specimen:

4 mL Lavender EDTA Vacutainer tube

Collection Information:

A minimum volume of 3 mL is required for testing.

Do not freeze the specimen. Send specimen in original collection tube. Do not transfer to polypropylene transfer vial.

Note: This test is only available to gastroenterologists at SJH/LHSC. If requested by other physician's approval must be given by the Biochemist on-call.

Do not transfer to polypropylene transfer vial.

Reference Ranges:

Results interpreted on report

Storage and Shipment:

Please send on ice packs. Specimen must be received at the testing site within 7 days of collection. If your test order does not specify Genotyping (or use the correct test mnemonic) TPMT Activity will be analyzed.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)