

TOXOPLASMA
SEROLOGY -
TOXOPLASMA
IGG/IGM

Orderable - TOXOG/TOXOM

Turn Around Time: 48 hours

Specimen:

5 mL Gold or 6 mL Red top Vacutainer tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

See report
Indeterminate Toxo. IgG results will be forwarded to Public Health for confirmation on weekdays.

Interpretive Comments:

A single serum sample will yield significant results for establishing previous exposure only.

Critical Information Required:

Provide clinical history.

Indicate purpose of the request as diagnosis of acute disease (IgM) or determining previous exposure (IgG).



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Electrochemiluminescence
Immunoassay



Test Schedule:
24/7