

TRAZODONE

Orderable - TRAZ

Turn Around Time: 10 days

Alternate Name(s):

Desyrel

Specimen:

6 mL Dark Green (Lithium Heparin) top Vacutainer tube

Collection Information:

Avoid gel-separator tubes

Minimum volume required is 1 mL

To monitor therapy, collect trough level specimen prior to morning dose or 10 - 12 h after last drug administration. Separate as soon as possible.

Reference Ranges:

2.2 - 4.3 $\mu\text{mol/L}$

Storage and Shipment:

Separate serum as soon as possible.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)