

TRIMIPRAMINE
SERUM/
PLASMA

Orderable - TRIM

Turn Around Time: 7 days

Alternate Name(s):

Surmontil

Specimen:

2 x 6 mL Red top Vacutainer tube
or
2 x Lavender top Vacutainer tubes

Collection Information:

Avoid gel-separator tubes

Minimum volume of serum or EDTA plasma required is 3 mL.

To monitor therapy, draw trough specimen 1 hour prior to morning dose or 10 -12 hours after last drug administration. Separate and freeze as soon as possible.

Reference Ranges:

170-1000 nmol/L

Toxic: Greater than 1700 nmol/L

Storage and Shipment:

Send frozen sample ASAP.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
GC/NP Detector



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)