

PNEUMOCOCCAL
ANTIBODY
TITRE

Orderable - PNEUMA(Pre), PNEUMB (post), PNEUM

Alternate Name(s):

Pneumococcal Antibody Titre, Pre Vaccine
Pneumococcal Antibody Titre, Post Vaccine
Pneumococcal Antibody Titre, Single Sample

Specimen:

Adult	Pediatric
6 mL Red top Vacutainer tube	0-2 years: 6 x 0.5 mL Red top Microtainer 2-10 years: 2 x 3 mL Red top tube

Collection Information:

3 mL serum required. Must indicate "pre" or "post".
Serum from gold top tube not acceptable.

Note: This test is only available at the request of Dr. Ashley Geerlinks, Dr. Kara Robertson Dr. April Price, Dr. Mark Kuprowski, Dr. Harold Kim, Dr. Rongbo Zhu and Dr. Samira Jeimy.

Reference Ranges:

Results sent directly to physician. See report.

Special Processing:

Aliquot serum and freeze. Must send both Pre and Post specimens at the same time.

Comments:



Laboratory:
Core Lab



Requisition:
[SICK KIDS](#)
[PNEUMOCOCCAL](#)
[ANTIBODY TITRE](#)
[REQUISITION](#)



Test Schedule:
Referred out as required.
Pre and Post specimens
must be sent together.



Referred Out Location:
The Hospital for Sick
Children
Rapid Response
Laboratory- Room 3642
555 University Avenue,
Toronto, ON
M5G 1X8
Canada
Phone: 416-813-7200
Toll Free: 1-855-381-3212



PNEUMOCOCCAL
ANTIBODY
TITRE

Note: This test is only available through consultation with a Clinical Immunologist (Drs. Kim/Jeimy/Kuprowski) for the evaluation of primary immunodeficiency.

Storage and Shipment:

Freeze serum at -20°C. Ship on dry ice.