



ALPHA-
GALACTOSIDASE,
LEUKOCYTE

Orderable - AGALL

Alternate Name(s):

Fabry's Disease

Specimen:

6 mL Dark Green (Sodium Heparinized) top Vacutainer tube

Collection Information:

Sample must be collected so that it will be received in lab no later than 3:00 pm of collection day.

The preferred specimen is blood. Collect whole blood and send immediately. Sample must be kept cold but not frozen.

Reference Ranges:

Leukocyte:
24-56 nmol/hr/mg protein

Storage and Shipment:

Send to lab **immediately**.



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
As required