

ALPHA-
GALACTOSIDASE,
PLASMA

Orderable - AGALP

Alternate Name(s):

Fabry's Disease

Specimen:

6 mL Dark Green (Sodium Heparinized) top Vacutainer tube

Collection Information:

Minimum volume required is 1 mL plasma sent frozen.

Reference Ranges:

4.2-15 nmol/hr/mL plasma

Special Processing:

Spin and freeze plasma -20°C.



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
As required