

ALPHA-N-
ACETYLGALACTOSAMINIDASE,
LEUKOCYTE
/PLASMA/
FIBROBLASTS

Orderable - GALNACL/GALNACP/ Fibroblast - Use requisition to order

Alternate Name(s):

Schindler Disease

Specimen:

1. 2 x 6 mL Dark Green (Sodium Heparinized) top Vacutainer
2. 4.5 mL Green (Lithium Heparin) top Vacutainer tube (preferred specimen type)
3. Fibroblasts

Collection Information:

Sample must be collected so that it will be received in lab no later than 3:00pm of collection day.

1. Send 8-10 mL whole blood immediately. Keep cold but not frozen.
2. The preferred specimen is plasma. A minimum volume of 1 mL of plasma is required for testing.
3. Collect fibroblasts and send in a sterile container.

Reference Ranges:

Leukocyte:
0.55-1.10 deltaA410/hr/mg protein

Fibroblast:
1.83-2.67 deltaA410/hr/mg protein

Storage and Shipment:

Sample must be sent to lab **immediately**.



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
As required