

# AMINO ACIDS, 24-HOUR URINE

## Orderable - AAU24

Turn Around Time: 3 days

### Alternate Name(s):

Cystine  
Cystinuria

### Specimen:

24-hour urine collection

### Collection Information:

#### Collection of 24-Hour Urine Sample

No Preservative necessary, but specimen should be kept refrigerated during collection.

Aliquot 10 mL. Store at -20°C and send frozen to Biochemical Genetics.

Please complete the laboratory label attached to the collection bottle to include patient's name, collection start time and date and collection end time and date.

### Reference Ranges:

See report

### Storage and Shipment:

Keep specimen refrigerated during collection. Freeze a 10 mL aliquot at -20°C and send frozen to Biochemical Genetics.



**Laboratory:**  
Biochemical Genetics Lab



**Requisition:**  
[BIOCHEMICAL GENETICS  
LAB REQUISITION](#)



**Test Schedule:**  
As required