

AMINO ACIDS, 24-HOUR URINE

Orderable - AAU24

Turnaround Time: 5 days

Alternate Name(s):

Cystine
Cystinuria

Specimen:

24-hour urine collection

Collection Information:

[Collection of 24-Hour Urine Sample](#)

No Preservative necessary, but specimen should be kept refrigerated during collection.

Aliquot 10 mL. Store at -20°C and send frozen to Biochemical Genetics.

Please complete the laboratory label attached to the collection bottle to include patient's name, collection start time and date and collection end time and date.

Reference Ranges:

See report

Storage and Shipment:

Keep specimen refrigerated during collection. Freeze a 10 mL aliquot at -20°C and send frozen to Biochemical Genetics.



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
As required