

AMINO ACIDS, RANDOM URINE

Orderable - AAU

Turn Around Time: 3 days

Specimen:

Random urine

Collection Information:

Volume of urine required is 2 mL.

Store urine at -20°C and send frozen to Biochemical Genetics.

Reference Ranges:

See report

Special Processing:

Send urine sample frozen.

Storage and Shipment:

Store urine at -20°C and send frozen to Biochemical Genetics.



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
As required