

GALACTOSE-1-  
PHOSPHATE  
URIDYLTRANSFERASE -  
QUALITATIVE,  
DRIED BLOOD  
SPOT

## Orderable - G1PUDBS

### Alternate Name(s):

Galactosemia Screen

### Specimen:

Dried blood spots on a Newborn Screening-type card.

Blood can be from a heel prick, finger prick, or Green top (Sodium or Lithium Heparinized) Vacutainer tube.

### Collection Information:

On a Newborn Screening-type card, saturate each circle completely (i.e blood dot should be clearly visible on both front and back of card).

Air dry horizontally at room temperature. Do **NOT** heat the card to dry!

Once the blood spots are completely dry (approximately 3 hours) seal the card in a plastic bag and send to Biochemical Genetics at room temperature.

### Reference Ranges:

See report



**Laboratory:**  
Biochemical Genetics Lab



**Requisition:**  
[BIOCHEMICAL GENETICS  
LAB REQUISITION](#)



**Test Schedule:**  
As required