

METHYLMALONIC  
ACID,  
URINE

## Orderable - MEMAU

Turn Around Time: 30 days

### Alternate Name(s):

Methylmalonic Aciduria Monitoring

### Specimen:

Random urine

### Collection Information:

Collect a minimum volume of 10 mL of urine (random) in a sterile container, then freeze.

Keep at -20°C

### Reference Ranges:

See report

### Storage and Shipment:

Freeze sample and send to lab **immediately**.



**Laboratory:**  
Biochemical Genetics Lab



**Requisition:**  
[BIOCHEMICAL GENETICS  
LAB REQUISITION](#)



**Test Schedule:**  
As required