

METHYLMALONIC ACID, URINE

Orderable - MEMAU

Turnaround Time: 40 days

Alternate Name(s):

Methylmalonic Aciduria Monitoring

Specimen:

Random urine

Collection Information:

Indicate date of birth and clinical diagnosis. This test may not be done if a preliminary Organic Acid test is normal.

Collect a minimum volume of 10 mL of urine (random) in a sterile container, then freeze.

Reference Ranges:

Qualitative Test: Normal

Storage and Shipment:

Store and send frozen.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Gas Chromatography with
Mass Spec.



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)