

ORGANIC ACID – SCREENING, URINE

Orderable - ORGA

Note:

Methylmalonic acid urine is part of urine organic acids testing.

Specimen:

Random urine

Collection Information:

Collect a minimum volume of 5 mL of urine (random) and send frozen to lab.

Keep at -20°C

Reference Ranges:

See report

Storage and Shipment:

Send frozen sample to lab **immediately**.



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
As required