

OROTIC ACID, URINE

Orderable - OROTO

Alternate Name(s):

Urea Cycle Disorder

Specimen:

Random urine

Collection Information:

Collect a minimum volume of 10 mL of urine (random) in a sterile container and send frozen to lab.

Keep at -20°C.

Reference Ranges:

0.0-7.7 ug/mg creatinine

Storage and Shipment:

Send frozen sample to lab **immediately**.



Laboratory:
Biochemical Genetics Lab



Requisition:
[BIOCHEMICAL GENETICS
LAB REQUISITION](#)



Test Schedule:
As required