

CHRONIC
GRANULOMATUS
SCREENING
INVESTIGATION

Orderable - FLOW

Turnaround Time: 24 hours STAT: 4 hours

Alternate Name(s):

Oxidative Burst

Specimen:

Adult	Pediatric
4 mL K2 or K3 EDTA Lavender top Vacutainer tube	0-2 years: 0.5 mL Lavender Microtainer 2-10 years: 3 mL Lavender top

Collection Information:

Minimum volume required for testing is 0.5 mL.

Notify the laboratory (519-685-8300 ext. 57450) prior to sending sample.

For referred out testing please contact the Laboratory or Dr. Ben Hedley 519 685 8300 ext. 55402.

Reference Ranges:

Neutrophils demonstrate normal oxidative burst following stimulation by PMA, not suggestive of a diagnosis of CGD.

Interpretive Comments:

Negative: Neutrophils demonstrate normal oxidative burst following stimulation by PMA, not suggestive of a diagnosis of CGD

Positive (Afflicted): Neutrophils demonstrate abnormal oxidative burst following stimulation by PMA, suggestive of a diagnosis of CGD. CGD is an extremely rare disease. It is recommended that the test be repeated for confirmation with a new specimen.

Positive (Carrier): Two populations of neutrophils identified. One population demonstrates normal oxidative burst following stimulation by PMA, the second demonstrates abnormal oxidative burst. Suggests patient is a carrier for CGD. CGD is



Laboratory:
Flow Cytometry Lab



Requisition:
GENERAL LABORATORY
REQUISITION




Method of Analysis:
Flow cytometry



Test Schedule:
Monday to Thursday only
0800-1300.
Do not send samples on
Fridays.

Prior to STAT holidays,
routine specimens must
be received in the Flow
Cytometry Laboratory at
Victoria Hospital by 1300.



CHRONIC
GRANULOMATOUS
SCREENING
INVESTIGATION

an extremely rare disease. It is recommended that the test be repeated for confirmation with a new specimen.

Special Processing:

Prior booking of the testing is required. Call 519-685-8300 x 57450

Storage and Shipment:

Ship sample at room temperature (20°C to 25°C) prior to testing.

Samples are to be shipped immediately after collection.

Samples must be received for testing within 4 hours of collection.

Specimen Stability:

Sample is stable for 6 hours at room temperature.