

SCREEN FOR
FACTOR V
LEIDEN
MUTATION

Orderable - FVL

Turnaround Time: 2 weeks

Alternate Name(s):

FVL

Specimen:

Adult	Pediatric
4 mL Lavender (EDTA) top Vacutainer tubes	0-2 years: Contact Lab at Ext. 58140 for number of tubes required prior to sampling

Collection Information:

LHSC patient samples require authorization by a Hematologist.

Please send the specimen unspun and at room temperature.

Reference Ranges:

Negative



Laboratory:
Molecular Diagnostics Lab



Requisition:
[MOLECULAR DIAGNOSTIC
REQUISITION](#)



Method of Analysis:
Automated clot detection



Test Schedule:
Bi-weekly,
Monday to Friday 0800-
1600 hours