

SCREEN FOR  
FACTOR V  
LEIDEN  
MUTATION

## Orderable - FVL

Turnaround Time: 2 weeks

### Alternate Name(s):

FVL

### Specimen:

Adult	Pediatric
4 mL Lavender (EDTA) top Vacutainer tubes	0-2 years: Contact Lab at Ext. 58140 for number of tubes required prior to sampling

### Collection Information:

LHSC patient samples require authorization by a Hematologist.

Please send the specimen unspun and at room temperature.

### Reference Ranges:

Negative



**Laboratory:**  
Molecular Diagnostics Lab



**Requisition:**  
[MOLECULAR  
DIAGNOSTICS  
REQUISITION](#)



**Method of Analysis:**  
Automated clot detection



**Test Schedule:**  
Bi-weekly,  
Monday to Friday 0800-  
1600 hours