

CLOZAPINE,  
SERUM/  
PLASMA

## Orderable - CLOZ

Turnaround Time: 5 days

### Alternate Name(s):

Clozaril

### Specimen:

4 mL K<sub>2</sub> or K<sub>3</sub> EDTA Lavender top Vacutainer tube

### Collection Information:

Serum from a 6 mL Red top Vacutainer tube is also acceptable.

No Gel Separator Tubes

Minimum volume required for testing is 0.5 mL.

### Reference Ranges:

Minimum effective concentration is 1070 nmol/L

### Interpretive Comments:

Norclozapine, metabolite: not active.

Alert value is 3100 nmol/L

### Storage and Shipment:

Spin immediately, aliquot plasma or serum and freeze.

Ship on ice packs or frozen.



**Laboratory:**  
Toxicology/Special  
Chemistry Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
LC-MS/MS



**Test Schedule:**  
Monday – Friday  
0800-1600