

HALOPERIDOL

## Orderable - HALO

Turnaround Time: 5 days

### Alternate Name(s):

Haldol

### Specimen:

6 mL Red top or 4.5 mL Green (Lithium Heparin) top Vacutainer tube

### Collection Information:

Minimum volume required is 0.5 mL of serum or plasma.

Collect blood aseptically in a Vacutainer tube.

### Reference Ranges:

11-53 nmol/L

### Storage and Shipment:

Separate and freeze. Ship on ice packs or frozen.



**Laboratory:**  
Toxicology/Special  
Chemistry Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
LC-MS/MS



**Test Schedule:**  
Monday--Friday  
0800-1600