

METHADONE
URINE
QUALITATIVE

Orderable - METHAU

Turnaround Time: 5 days

Alternate Name(s):

Metadol

Specimen:

Urine

Collection Information:

Minimum 10 mL random urine collected in a sterile container.

Reference Ranges:

See report

Storage and Shipment:

Keep sample cold or frozen for shipment to laboratory.



Laboratory:
Toxicology/Special
Chemistry Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
GC-MS



Test Schedule:
Monday--Friday
0800-1600