





Pathology and Laboratory Medicine

Orderable - FVAL

Turnaround Time: 5 days

Alternate Name(s):

Depakene



Laboratory:

Toxicology/Special Chemistry Lab



Requisition:

GENERAL LABORATORY REQUISITION



Method of Analysis:

LC-MS/MS

Specimen:

 $2 \times 5 \text{ mL}$ Gold top Vacutainer tube

or

2 x 6 mL Red top Vacutainer tube

or

2 x 6 mL Dark Green (Sodium Heparin) top Vacutainer tube

Collection Information:

Minimum volume of serum required for testing is 2 mL.

Spin and aliquot serum or plasma and freeze immediately and send to the lab.

Reference Ranges:



Test Schedule: Monday--Friday

0800-1600

35.0-140.0 µmol/L

Storage and Shipment:

Separate and freeze. Ship on ice packs or frozen.