

MYASTHENIA  
GRAVIS  
PANEL  
(anti-AChR  
and MuSK)

## Orderable - ACHR/MUSK

Turnaround Time: 7 days

### Alternate Name(s):

ACHR Antibodies  
AChRab  
MuSK Antibodies  
MuSK Ab



**Laboratory:**  
Clinical Immunology Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Cell-based assay (CBA)



**Test Schedule:**  
Batch analysis

### Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube or 6 mL Red top Vacutainer tube	0-2 years: 0.5 mL Red top Microtainer 2-10 years: 2 mL Red top Vacutainer tube
EDTA, heparin or citrated plasma are also acceptable	

### Collection Information:

Minimum volume of serum is 1 mL (2 mL preferred)

### Reference Ranges:

Negative

### Interpretive Comments:

Antibodies associated with Myasthenia Gravis.

### Comments:

This test is restricted to neurologists and ophthalmologists at LHSC/SJHC. Accepted from referred in locations.

MYASTHENIA  
GRAVIS  
PANEL  
(anti-AChR  
and MuSK)

**Storage and Shipment:**

Sample can be stored for up to 14 days at 4°C. For longer storage, store at -20°C.